

Anatoxin-a, Cylindrospermopsin, Microcystin & Saxitoxin Report
(Utah Division of Water Quality)

<u>Sample Identification</u>	<u>Site</u>	<u>Collection Date</u>
4994790	Jordon R @ Utah Lake Outlet	10/22/2014
4994720	Jordon R @ Narrows Pump Station	10/22/2014
NA	Middle of Lindon Harbor	10/22/2014
NA	East Side of Lindon Harbor	10/22/2014
NA	Middle of Utah Lake State Park Harbor	10/22/2014

Toxins – Anatoxin-a (ANTX-A), cylindrospermopsin (CYN), microcystin (MC), saxitoxin (STX)

Sample Prep – The samples were ultra-sonicated to lyse all cells and release toxins. Solid phase extraction (SPE) was also utilized for anatoxin-a prep and pre-concentration (20x) followed by filtration. Lab Fortified Matrix (LFM) duplicates were prepared at 1.0 µg/L ANTX-A, 1.0 µg/L CYN and 0.2 µg/L STX.

Analytical Methodology – Liquid chromatography/ mass spectrometry/ mass spectrometry (LC/MS/MS) was utilized for the determination of ANTX-A. The $[M+H]^+$ ion for ANTX-A (m/z 166) was fragmented and the major product ions (m/z 91, 106, 131, 149) were monitored. The current method limit of detection (LOD) is 0.05 µg/L, with a limit of quantification (LOQ) of 0.1 µg/L for ANTX-A.

A microcystins enzyme linked immunosorbent assay (ELISA) was utilized for the quantitative and sensitive congener-independent detection of MCs. The current assay is sensitive to down to a LOD/LOQ of 0.15 µg/L for total MCs. The average recoveries of laboratory fortified blanks (LFB) spiked with 1 µg/L MCLR were 104%.

A cylindrospermopsin enzyme linked immunosorbent assay (ELISA) was also utilized for the quantitative detection of CYN. The current assay is sensitive down to a LOD/LOQ of 0.1 µg/L for CYN. A lab fortified blank (LFB) spiked with 1.0 µg/L CYN was recovered at 106%.

A saxitoxin enzyme linked immunosorbent assay (ELISA) was utilized for the quantitative detection of saxitoxin. The current assay is sensitive down to a LOD/LOQ of 0.02 µg/L saxitoxin. The LFB (0.2 µg/L STX spike) recovery was 85%.

Summary of Results

(µg/L)

<u>Sample</u>	<u>MC</u> (ELISA)	<u>CYN</u> (ELISA)	<u>STX</u> (ELISA)	<u>ANTX-A</u> (LC-MS/MS)
Jordon R @ Utah Lake Outlet	ND	ND	ND	ND
Jordon R @ Narrows Pump Station	0.19	ND	ND	ND
Middle of Lindon Harbor	ND	ND	ND	ND
East Side of Lindon Harbor	3.05	ND	ND	ND
Middle of Utah Lake State Park Harbor	ND	ND	ND	ND
<i>Detection Limits (µg/L)</i>	<i>0.15</i>	<i>0.10</i>	<i>0.05</i>	<i>0.05</i>

ND = Not detected above the detection limit

Submitted by:



Mark T. Aubel, Ph.D.

Date:

10/24/14